



38-50 Bedford Street
GILLMAN SA 5013
AUSTRALIA
Phone: (61 8) 8440 2000
Fax: (61 8) 8440 2001
www.chemicalreagents.com.au

SPECIFICATION SHEET

PRODUCT	<i>Australian Chemical Reagents</i> Standard Copper Solution 1000 mg/L Cu (Cyanide Matrix)
CATALOGUE NO	3353

PRODUCT SPECIFICATION:

Description: Standard Copper Solution
Nominal Concentration: 1000 ug / ml (1000 ppm)
Starting Material: Copper shot 99.999%
Matrix: alkaline cyanide 0.5% NaCN
Accuracy: $\pm 0.5\%$ of nominal concentration

Australian Chemical Reagents (ACR) does not warrant that this product is suitable for any use or purpose. The user must ascertain the suitability of the product for any intended purpose. Preliminary testing of the product before use or application is recommended. Any reliance or purported reliance upon ACR with respect to any skill or judgment or advice in relation to the suitability of this product for any purpose is disclaimed. Except to the extent prohibited at law, any condition implied by any statute as to the merchantable quality of this product or fitness for any purpose is hereby excluded. This product is not sold by description. Where the provisions of Part V, Division 2 of the Trade Practices Act apply, the liability of ACR is limited to the replacement or supply of equivalent goods or payment of the cost of replacing the goods or acquiring equivalent goods.