

PRODUCT SPECIFICATION SHEET

Spec. ACID - 01/07/01

Product Name:	Decon ACID RINSE
Description:	Acid rinse (neutraliser) and/or preliminary cleaner (pre-wash), for use in laboratory, medical and industrial automatic glassware washing machines. NOT suitable for <i>domestic</i> applications. Supplied as a liquid concentrate, for use in automatic machines.
Properties:	Biodegradable*, bactericidal**, non-flammable, non-foaming and totally rinsable.
Appearance:	Colourless, non-viscous liquid, having mild odour.
pH:	Decon ACID RINSE is acidic. pH 3.3 - 1% solution in water (typical).
Packaging:	Supplied in <i>5 litre</i>, high density polyethylene containers, with tamper-evident screw-on caps.
Applications:	For use in automatic washing machines as an acid rinse (neutraliser) and/or a pre-wash (preliminary cleaner), when cleaning/decontaminating glassware, ceramics, rubbers, plastics, stainless steel and ferrous metals. Being acidic, Decon ACID RINSE is NOT suitable for use on <i>non-ferrous metals</i>, notably <i>aluminium</i> and <i>zinc</i>.
Formulation:	A blend of aliphatic polycarboxylic acid, aliphatic alcohol and non-ionic surface active agent in an aqueous base.

*** Ready Biodegradability by OECD 301E (ISO method 7827 - 1986 (E))**

**** Satisfies the criteria in BS EN 1040 - 1997 for bacterial activity.**