

Dimethylformamide, Anhydrous (100 ppm)

HCON(CH ₃) ₂	FW. 73.10	Density 1 L = 0.949 Kg.
CAS-No.	68-12-2	Melting Point - 61 °C
Code	AH1053	Boiling Point 153 °C

Specifications

Assay (by GC.)	99.9%	min.
Identity (IR)	Passes test	
Color (APHA)	10	max.
Water (by Coulometry)	0.01%	max.
Acidity (mEq./g.)	0.0005	max.
Residue on Evaporation	0.001%	max.

Added with Molecular sieve.