

PRODUCT SPECIFICATION SHEET

Spec. DMATIC - 01/07/01

| | |
|----------------------|--|
| Product Name: | Decon DECOMATIC |
| Description: | Surface active cleaning agent, for use in laboratory, medical and industrial automatic glassware washing machines. NOT suitable for <i>domestic</i> applications. Supplied as a liquid concentrate, for use in automatic machines. |
| Properties: | Biodegradable*, bactericidal** phosphate-free, non-foaming, non-flammable and totally rinsable. |
| Appearance: | Clear to pale amber, non-viscous liquid, having faint odour. |
| pH: | Decon DECOMATIC is <i>alkaline</i>. The concentrate has a pH in excess of pH 13. |
| Packaging: | Supplied in <i>5 litre, 10 litre and 20 litre</i> high density polyethylene containers, with tamper-evident screw-on caps. |
| Applications: | For machine cleaning and/or decontaminating glassware, ceramics, rubbers, most plastics (NOT <i>polycarbonate</i>), stainless steel and ferrous metals. Being alkaline, Decon DECOMATIC is NOT suitable for use on <i>non-ferrous metals</i>, notably <i>aluminium</i> and <i>zinc</i>. |
| Formulation: | A blend of amphoteric surface active agents, non-phosphate detergent builders, alkalis and sequestering agents, in an aqueous base. |

* Ready Biodegradability by OECD 301E (ISO method 7827 - 1986 (E))

** Satisfies the criteria in BS EN 1040 - 1997 for bacterial activity.